

Corning® hybrigro SF™ Medium for High Density Hybridoma Culture and Increased Production

CORNING

SnAPPShots

A brief technical report from the Corning Applications Group

Hilary Sherman and Mark E. Rothenberg Ph.D. Corning Incorporated, Life Sciences Kennebunk, Maine

Introduction

There is potential that the use of antibodies as therapeutic candidates will increase significantly over the next few years. With this increase, more efficient and cost effective methods for manufacturing antibodies will become necessary. Corning® hybrigro SF™ medium aims to solve this problem by increasing hybridoma yields and production efficiency without the use of costly serum. In this paper we demonstrate how hybrigro SF outperforms two commercially available serum-free hybridoma media, as well as a commonly used serum-containing formulation resulting in higher cell density and protein production using two different hybridoma cell lines.

Method and Materials

MH677 (proprietary hybridoma cell line) and AE1 (ATCC Cat. No. HB-72) cells were cultured in Dulbecco's Modification of Eagle's medium (DMEM) (Corning Cat. No. 10-013-CM) supplemented with 10% fetal bovine serum (FBS) (Corning Cat. No. 35-010-CV)

and enumerated using the BioProfile® Flex analyzer (Nova Biomedical). MH677 cells were resuspended in either hybrigro SF medium (Corning Cat. No. 40-215-CV) or competitor serum-free hybridoma media. AE1 cells were resuspended in either hybrigro SF medium, DMEM supplemented with 10% FBS, or competitor serum-free hybridoma medium. Cells were seeded into 125 mL Erlenmeyer flasks (Corning Cat. No. 431405) at a concentration of 5×10^4 cells/mL. Cells were grown for 4 days in a humidified incubator at 5% CO₂ and 37°C at 90 rpm. Daily culture samples were assessed for cell viability, nutrient, and metabolite analysis on the BioProfile Flex. Samples were also collected from each vessel for protein production (Mouse IgG2a), which was assessed by following the IgG2a ELISA protocol provided by Alpha Diagnostic International (Alpha Diagnostic International Cat. No. 6340) for MH677 or the Easy Titer IgG Assay kit from Pierce (ThermoFisher Cat. No. 23300).

Results

MH677 Cell Growth and Function

MH677 cells cultured in the hybrigro SF medium exhibited higher cell density than cultured in either of the two serum-free hybridoma media tested. The analysis shows the results on the fourth day of growth to be statistically significant (ANOVA Newman-Keuls Post Test) (Figure 1). Additionally, MH677 cells grown in hybrigro SF medium had statistically higher IgG2a production as compared to both competitor serum-free hybridoma media (ANOVA Newman-Keuls Post Test) (Figure 2).




Figure 1. MH677 cells showed significantly increased cell densities by day four compared to competitor serum-free hybridoma media (n=9 ANOVA Newman-Keuls Post Test).




Figure 2. MH677 cultured in Corning hybrigro SF showed statistically higher levels of IgG2a production compared to both competitor serum-free media (n=9 ANOVA Newman-Keuls Post Test).




Figure 3. AE1 cells cultured in Corning hybrigo SF medium showed significantly increased cell densities by day four compared to competitor serum-free medium and a serum containing medium (n=9 ANOVA Bonferroni Post Test).

AEI Cell Growth and Function

AEI cells grown in Corning® hybrigo SF™ medium not only outperformed cells grown in a commercially available serum-free hybridoma medium, but they also outperformed cells growing in serum-containing medium in terms of cell density achieved after 4 days of growth (Figure 3). This was found to be statistically significant (ANOVA).

In addition to cell growth, AEI cells cultured in hybrigo SF medium also had statistically significant higher total protein production after 4 days of growth (ANOVA) compared to serum-containing and competitor serum-free media (Figure 4).

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, please visit www.corning.com/lifesciences/media or call 1.800.492.1110. Outside the United States, please call 978.442.2200.

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

[www.corning.com/
lifesciences/media](http://www.corning.com/lifesciences/media)

Worldwide Support Offices

ASIA/PACIFIC
Australia/New Zealand
t 0402-794-347

China
t 86 21 2215 2888
f 86 21 6215 2988

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries

t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA

Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589




Figure 4. AE1 cells cultured in Corning hybrigo SF medium produced significantly higher protein production after four days of culture compared to a competitor serum-free medium and a serum-containing medium (n=9 ANOVA Bonferroni Post Test).

Conclusions

- ▶ Higher MH677 and AE1 cell densities can be achieved with hybrigo SF medium as compared to other commercially available serum-free hybridoma media.
- ▶ MH677 cells and AE1 cells produced more mouse IgG when grown in hybrigo SF medium as compared to equivalent serum-free hybridoma media.
- ▶ AE1 cells achieved higher densities and higher protein production in hybrigo SF medium as compared to DMEM containing 10% FBS.
- ▶ hybrigo SF medium is an ideal choice for proliferation of hybridoma cultures and the production of monoclonal antibodies.