

Certificate of Analysis

Description: McCoy's 5A, 1X

(Iwakata & Grace Mod.)
with L-glutamine

Catalog #: 10-050-CV

Lot #: 10050091

Storage: Store refrigerated at 2-8°C

Expiration Date: 01/14

Country of Origin: USA

Notes:

Mediatech, Inc. products are prepared by a validated aseptic sterile filtration process to ensure that all products comply with the industry's Sterility Assurance Level (SAL) of at least 10^{-3} (i.e. demonstrating that a product produced in this manner has no more than 1 random contaminant per 1000 units).

This product was tested in accordance with currently approved Mediatech, Inc. specifications and procedures. Testing procedures are maintained in compliance with the current versions of the USP and/or EP, where applicable.

For IVD, Manufacturing & Research Uses Only. Utilization of this product apart from the labeled intended use may be a violation of local and/or Federal Law.

TEST PARAMETERS	SPECIFICATION	RESULT
pH	7.2 ± 0.2	7.1
Osmolality	290 ± 30 mOsm/Kg H ₂ O	297 mOsm/Kg H ₂ O
Endotoxin	≤ 0.25 EU/mL @FW	< 0.0010 EU/mL
Mycoplasma	Tested Negative	Negative
Sterility	Pass	Pass
Cell Line		
L929 (Cytotoxicity)	Pass	Pass
LM (Growth Promotion)	Pass	Pass
PRK (Growth Promotion)	Pass	Pass
NG-108 (Growth Promotion)	Pass	Pass

For detailed information on Mediatech, Inc. testing procedures, please contact our Technical Services department at (800) CELLGRO (235-5476) or via email at techserv@cellgro.com.

Written By/Date: *S. L. N. 01 Aug 2012*

Verified By/ Date: *[Signature] 01 Aug 2012*

1.800.CELLGRO TOLL FREE
703.471.5955 LOCAL AREA
703.471.1086 FACSIMILE

9345 DISCOVERY BOULEVARD MANASSAS, VIRGINIA