

Change Notification

Date:	22-DEC-2016
Material Description:	100 Liter Drum (Non-nestable)
Material Number:	M 05201
Change Description:	Drum Change


Dear Valued Customer:

In order to serve you better, Mediatech is updating the drum container used to store and ship 100 L bag products. The following table outlines the changes:

	Current	New
Total Weight (pounds)	14.6	16.5
Wall Thickness (inches)	0.125	0.150
Container Height (inches)	29.1	29.375
OD (inches)	19.3	18.5
Opening ID (inches)	17.5	16.87

The new drum is constructed of the same material and meets the same regulatory standard, UN 1H2/Y180/S**, as the current drum. In addition, the new drum has a flared base making it less prone to tipping. Please see Appendix A for products in scope for this update. The new drum will go into effect after current drum supplies are depleted. Manufactured lots will not use commingled new/old drum types.

For further information, contact Corning Scientific Support at 800-492-1110 or ScientificSupport@corning.com.


23 JAN 2017

Jason Jens
Plant Quality Manager
p (703) 736-1406
jensj@corning.com


23 Jan 2017

Andria Collier
Product Line Manager
p (703) 736-1428
collierap@corning.com



Mediatech, Inc.
A Corning Subsidiary
9345 Discovery Boulevard
Manassas, VA 20109
1-800-492-1110
www.cellgro.com

Appendix A

	Item ID	Description
1	25-055-LG	100L, Cell Culture Grade Water
2	25-065-LG	100L, Cell Culture Grade Water
3	25-065-LGR	100L, Cell Culture Grade Water