

Change Notification

| | |
|------------------------------|---------------------------------|
| Date: | 10-FEB-2017 |
| Material Description: | Molecular Biology Grade Water |
| Catalog Number(s): | 46-000-CI, 46-000-CM, 46-000-CV |
| Change Description: | Specification Update |


Dear Valued Customer:

In an effort to continuously improve the product provided to you, Mediatech is updating the pH specification for products with the prefix 46-000, Molecular Biology Grade Water. The following table outlines the change:

| | Current | New |
|-----------|--------------------|----------------------|
| pH | 6.0 ±1.0 @ 20-25°C | 6.0 ±1.0 @ 25°C ±1°C |

This update will align our pH specification with the current United States Pharmacopoeia (USP) Chapter <645>, *Water Conductivity*, Bulk Water Stage 3 testing acceptance criteria. This change better serves our customers by aligning with scientifically-validated USP requirements.

The update is effective moving forward and impacts products with the SKU prefix 46-000. Current supplies of 46-000-XX will retain the current specification. There has been no change to the product manufacturing process. For further information, contact Corning Scientific Support at 800-492-1110 or ScientificSupport@corning.com.



13 Feb 2017

Jason Jens
Plant Quality Manager
p (703) 736-1406
jensj@corning.com



21 Feb 2017

Andria Collier
Product Line Manager
p (703) 736-1428
collierap@corning.com



Mediatech, Inc.
A Corning Subsidiary
9345 Discovery Boulevard
Manassas, VA 20109
1-800-492-1110
www.cellgro.com