

CHANGE NOTIFICATION LETTER
OSMOLALITY SPECIFICATION CHANGE FOR 10-027
DECEMBER, 2015

The osmolality specification for catalog number 10-027 will change from 315 ± 20 mOsm to 305 ± 15 mOsm.

This change in osmolality specification is based on historical data and is not associated with any product formulation changes. The revised specifications are tighter and correlate better to the mathematically derived osmolality calculations based on formula.

The change will go into effect immediately and all newly manufactured lots of 10-027 will be required to meet the revised specifications.

Mediatech and Corning are committed to providing our customers with consistently high quality products, and we are making these changes in support of this commitment.

If you have any questions regarding this change, please contact Corning Scientific Support (Scientificsupport@corning.com or (800) CELLGRO (235-5476)).