

CHANGE NOTIFICATION LETTER

Foam Change for 100L Drums

January, 2016

You are receiving this notification because a product custom manufactured for you is affected by this change.

The padding for the 100L drums will change from a polyurethane (PU) material to a polyethylene (PE) material. Mediatech has used PE foam void filler in other products with no negative effect. The PE foam pads will provide the same void filler and protection as the current PU foam pads and are the same thickness. The change from PU to PE foam has several advantages.

1. PE foam is cleaner than PU foam with a more open cell by design, allowing for less collection of debris and moisture within the foam.
2. The PE foam is less flammable than the PU foam.

The change is planned to take place for products made on or after 15Mar2016.

If you have any questions regarding this change, please contact Corning Scientific Support (Scientificsupport@corning.com or (800) CELLGRO (235-5476)).

The logo consists of the word "CORNING" in a white, sans-serif font, centered within a solid orange square.

Mediatech, Inc.
A Corning Subsidiary
9345 Discovery Boulevard
Manassas, VA 20109
Main (703) 471-5955
Customer Service (800) 235-5476
www.cellgro.com