

Corning® glutagro™ Medium and Supplements

CORNING

The Best Research Requires The Healthiest Cells.

Corning glutagro medium and supplement combine the quality you expect from Corning with the reagent stability you need. Corning glutagro medium and supplement are the perfect replacement for glutamine.

Cells grown in traditional media supplemented with L-glutamine are susceptible to cytotoxicity. With Corning glutagro medium and supplement, cells receive this needed supplement in the healthiest environment possible.

This modified form of L-glutamine does not degrade into ammonia during storage and incubation, minimizing cytotoxicity and improving culture quality.

Ordering Information


| Cat. No. | Description | Unit Size | Qty/Pk |
|-----------|--|-----------|--------|
| 25-015-CI | Corning glutagro medium and supplement | 100 mL | 1 |
| 10-101-CV | DMEM with glutagro supplement | 500 mL | 6 |
| 10-102-CV | DMEM without sodium pyruvate | 500 mL | 6 |
| 10-103-CV | DMEM/F-12 with glutagro supplement | 500 mL | 6 |
| 10-104-CV | RPMI 1640 with glutagro supplement | 500 mL | 6 |

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, please visit www.corning.com/lifesciences/media or call 1.800.235.5476.



www.corning.com/lifesciences/solutions

Mediatech, Inc.
A Corning Subsidiary
9345 Discovery Boulevard
Manassas, VA 20109
t 800.235.5476
t 703.471.5955
f 703.467.9851

www.corning.com/lifesciences/media

CORNING | **FALCON®** | **AXYGEN®** | **GOSSELIN™** | **PYREX®**

For a listing of trademarks, visit us at www.corning.com/lifesciences/trademarks. All other trademarks in this document are the property of their respective owners.